

Codice	PREPAZ05	
Titolo del documento	RACCOLTA DI FECI PER RICERCA SANGUE OCCULTO	
Area	Laboratorio	Rev. 2 del 10/11/2022

CONTENITORE

Le feci per la ricerca del sangue occulto devono essere raccolte attraverso il contenitore rilasciato dal Punto Prelievi, seguendo le annesse istruzioni. Il contenitore deve essere tenuto al di fuori della portata dei bambini. Evitare che il liquido entri in contatto con occhi e mucose.

RACCOMANDAZIONI

- Non è necessario praticare alcuna particolare dieta prima del test.
- Non eseguire il test durante il periodo mestruale, in caso di diarrea e in presenza di sanguinamento del retto (emorroidi)
- Evitare di contaminare le feci con le urine e l'acqua di scarico
- Si consiglia di disporre uno strato di carta igienica all'interno del W.C. in modo da trattenere le feci per il tempo necessario ad eseguire il prelievo.
- Se la prescrizione medica prevede la ricerca di sangue occulto fecale su tre campioni, i prelievi vanno effettuati in tre diverse sedute.


MODALITA' DI RACCOLTA

1. Identificare il contenitore con i propri dati anagrafici. Svitare il tappo verde con l'asta di prelievo, mantenendo il tubo inverticale. Non agitare il tubo capovolto.

2. Inserire l'asta nelle feci in quattro punti diversi e strisciare la croce sulla superficie. Ricoprire di feci solo le scanalature dell'asta. Non raccogliere feci in eccesso.

3. Reinserire il tappo verde e avvitare a fondo. Non riaprire. Agitare gentilmente il contenitore tramite inversione.

4. Restituire il dispositivo di raccolta per l'analisi entro una settimana a temperatura ambiente o mantenere il contenitore di raccolta in frigorifero e restituirlo entro 2 settimane.



Codice	PREPAZ05	
Titolo del documento	FECES COLLECTION FOR OCCULT BLOOD RESEARCH	
Area	Laboratorio	Rev. 2 del 10/11/2022

CONTAINER

Faeces for the search for occult blood must be collected through the container issued by the sampling point, following the attached instructions. The container must be kept out of the reach of children. Avoid liquid contact with eyes and mucous membranes.

RECOMMENDATIONS

- It is not necessary to practice any particular diet before the test.
- Do not perform the test during menstruation, in case of diarrhea and in the presence of bleeding from the rectum (haemorrhoids)
- Avoid contaminating feces with urine and wastewater
- It is advisable to place a layer of toilet paper inside the W.C. in order to hold the faeces for the time necessary to carry out the sampling.
- If the medical prescription provides for the search for faecal occult blood on three samples, the samples must be taken in three different sessions.

METHOD OF COLLECTION

- Identify the container with your personal data. Unscrew the green cap with the sampling rod, keeping the tube vertical.

Do not shake the tube upside down.

- Insert the rod into the stool at four different points and cross-swipe the surface.

Cover only the grooves of the shaft with faeces.
Do not collect excess feces.

- Reinsert the green cap and screw it tight

Don't reopen. Gently shake the container by inversion.

- Return collection device for testing within one week at room temperature or keep collection container refrigerated and return within 2 weeks.

